

How to Pair TeltoHeart with TeltoHeart mobile application

Start using TeltoHeart

How to Pair TeltoHeart with TeltoHeart mobile application



Tap Pair device to proceed with the pairing process.



TeltoHeart application will initiate a search for nearby devices.



Select your smart wearable by tapping on it.



TeltoHeart will notify you of an incoming pair request.



When prompted to do so, enter the passkey shown on your TeltoHeart into the Bluetooth pairing request window in the TeltoHeart application.

6

TeltoHeart will notify you about a successful pair with the TeltoHeart application.



7

Make sure TeltoHeart application is always running in the background, otherwise the data transfer will not be available.

Please note:



To ensure uninterrupted data transfer from your watch to the platform, it is important to keep the TeltoHeart application running in the background at all times. If you close the app, data transfer to the TeltoCare platform will not be available.